Investigation of PEALD SiO₂ thin films by rf-GDOES

Zhen Zhu^{1,2}, <u>Chiara Modanese</u>¹, Perttu Sippola¹, Marisa Di Sabatino³, Hele Savin¹

¹ Department of Electronics and Nanoengineering, Aalto University, Espoo, Finland ² Beneq Oy, Espoo, Finland

³ Dept. Materials Science and Engineering, Norwegian University of Science and Technology (NTNU), Trondheim, Norway



Aalto University School of Electrical Engineering



\Box NTNU

Norwegian University of Science and Technology

One Stop for All ALD





Development Services Coating Services

NAME AND

Industrial Equipment

Customer Services

Motivation - PEALD

- Thermal ALD (atomic layer deposition) is a thinfilm growth method with many advantages:
 - conformality
 - ultra-thin film
 - repeatability
 - ...
- PEALD (plasma enhanced ALD) provides further advantages:
 - It broadens applicable processing conditions (purge times, precursors, ...)
 - It opens for a wider range of material properties (thermally sensitive substrates, ...)
 - It shortens the deposition time







Aalto University School of Electrical Engineering

INTINU Norwegian University of Science and Technology



Motivation – PEALD-SiO₂

- SiO₂ thin-films are a well-studied material.
- Nevertheless, properties of SiO₂ deposited by PEALD under varying process parameters are not widely understood:
 - Effect of plasma power?
 - Effect of plasma temperature?
 - Effect of plasma composition?

alto University

ngineering

of Electrical

- Aim of our study:
 - tackle the challenge of determining the chemical composition of PEALD-SiO₂ thin-layer (d < 200 nm) to understand the effect of plasma power

Norwegian University of

Science and Technology

2. use a fast technique for routine analysis





TFS 200, Beneq

10/21/2019



Introduction – rf-GDOES

- rf-GDOES (radio frequency glow discharge optical emission spectroscopy) is an effective analytical technique for thin film analysis
- Several advantages:
 - very good depth resolution (~ nm)
 - wide dynamic range (12 orders) →
 detection limits (down to 1ppm)
 - multi-elemental analysis
 - fast analysis (few seconds-few minutes)
 - simple sample preparation
- Main limitation (not hindering analysis of ALD thin films):
 - mm-range lateral resolution



Photo and schematic diagram of the analytical GD source



Aalto University School of Electrical Engineering

NTNU Norwegian University of Science and Technology



Horiba, 2014

5

Materials



sputtered crater after GDOES measurement

- SiO₂ thin-films deposited on Czochralski Si substrates
- Deposition by O₂ based PEALD and under varying plasma powers:
 - 50 W
 - 180 W
 - 300 W



Aalto University School of Electrical Engineering Norwegian University of Science and Technology



Methods

- PEALD-SiO₂ deposition \rightarrow Beneq TFS 200 reactor ٠
- PEALD-SiO₂ compositional profiles (qualitative) \rightarrow rf-GDOES ٠ Horiba Profiler2:
 - plasma power: 35 W _
 - μs-pulsed discharge _
 - pre-cleaning and flushing (surface cleaning) _
 - total sputtering time < 20 s —



PEALD-SiO₂ thickness \rightarrow GDOES ? ٠

Electrical



www.horiba.com



Norwegian University of Science and Technology



10/21/2019

Qualitative elemental profiles

Example of elemental profiles for the PEALD-SiO₂ deposited at 180 W



→ The matrix and main impurity elements were successfully analysed



Aalto University School of Electrical Engineering



Chemical composition

Intensities of elemental emission lines were integrated over the SiO₂ thin-film depth



- Intensity difference in Si, O and H is within the accuracy of the technique (±10%)
- Clear intensity difference in C and N \rightarrow lower content with increasing plasma power ٠



ngineering

Aalto Universitv Electrical Norwegian University of Science and Technology



Structural properties

- Sputtering rate calculation:
 - same rate assumed for all films → sputtering time divided by the thickness measured by ellipsometry → calculation of the sputtering rate





Norwegian University of Science and Technology



10/21/2019

Structural properties

- Sputtering rate calculation:
 - same rate assumed for all films \rightarrow sputtering time divided by the thickness measured by ellipsometry \rightarrow calculation of the sputtering rate
- The sputtering rate decreases with increasing plasma power:



- film thickness measured by ellipsometry is similar for all films (148nm @50W, 146nm @180W and 145nm @300W) → reason for difference in sputtering time ?
- it seems to be influenced by film density \Leftrightarrow measured density increase of ~ 5% with increasing plasma power (\geq 180W) (from XRR)





Norwegian University of Science and Technology



• The thickness of a transparent thin film, *d* can be calculated as*:



$$d = \frac{k \times \lambda}{2 \times n}$$

- k is the number of oscillations of emission line i in the thin film
- λ is the wavelength of emission line *i*
- *n* is the refractive index of the thin film at wavelength *i*





Norwegian University of Science and Technology



 Comparison of thickness calculated from GDOES vs from reference techniques (ellipsometry and XRR):

Method	PE-ALD SiO ₂ thickness	PE-ALD plasma power		
		50W	180W	300W
GDOES	d _o (nm)	129	128	128
	d _{si} (nm)	129	128	127
Ellipsometry	d from ellipsometry (nm)	148	146	145
XRR	d from XRR (nm)	148	147	147

- → Note that $\frac{\lambda}{2 \times n}$ = 93 nm (corresponding to the thickness for k = 1) \Leftrightarrow good agreement with 100 nm value as suggested by Dorka *et al.**
- → Thus, for thin films (d < 200 nm), the relevance of the accuracy of k is dominant over the refractive index for the calculation of d.





NTNU Norwegian University of Science and Technology



*Dorka, Kunze, Hoffmann, JAAS 2000, 15, 873-876

Summary

- GDOES analysis provided important information on the chemical and structural properties of PEALD SiO₂ thin films, with good depth resolution and fast analysis.
- Detection limits are sufficient for the analytical purposes.
- GDOES analysis of thin (< 200 nm), transparent layers is challenging.
- The calculation of the thickness is severely affected by the number of oscillations k in the elemental lines, which cannot be rounded to the closest integer for small k values.







Thank you for your attention!

The authors acknowledge the European Research Council [EU's FP7 Programme ERC, Grant Agreement No. 307315] and the Micronova Nanofabrication Centre of Aalto University.



• The thickness of a transparent thin film, d can be calculated as*:

$$d = \frac{k \times \lambda}{2 \times n}$$

- *n* is the refractive index of the thin film at wavelength *i*

- Example of the simulated refractive index, n of the sample deposited at 180 W.
- The fitting is done on the measured values in the visible range.



• Influence of the number of oscillations within the film thickness:

SiO,	50W	180W	300W
k@130nm (O line)	4	4	4
n@130nm (O line)	2.018	2.035	2.030
k@288nm (Si line)	1.38	1.38	1.37
n@288nm (Si line)	1.545	1.552	1.556

Note that *n* is assumed to be constant throughouot the whole film thickness

→ Thus, for thin films (d < 200 nm), the relevance of the accuracy of k is dominant over that of the refractive index for the calculation of d.

Infrared spectrometry

- Qualitative analysis
- FTIR does not show clear differences among the thin films
- Transmission is also influenced by the Si substrate → limitation for possibility of evaluation of spectra



One Stop for All ALD





Development Services Coating Services

NAME AND

Industrial Equipment

Customer Services